


Steri-Vacc Transfer Set

Catalog Number: TS-917



Catalog Number: TS-917

Description: Hand operated 10 ml Syringe with 2-way stopcock pre-attached 4 mm PVC tubing for sterile docking.

Applications: For sterile docking to a primary bag of cells to be transferred to the AFC 6-F Steri-Vacc system. May be used for other sterile transfer applications.

Specifications/Technical Info:

One 10ml syringe with 2 way stopcock and 7 cm of 4 mm PVC tubing packaged individually sterile.

Features:

Easy to use.

Instructions for use.

Closed system compatibility.

Stopcock points to the "OFF" direction of flow.



American Fluoroseal Corporation.

431 -A East Diamond Ave. Gaithersburg, MD 20877
Tel: 301-990-1407 Fax 301-990-1472 www.toafc.com